BROOKS RAND LABS

MEANINGFUL METALS DATA

A Technique for the Simultaneous Determination of Monomethyl and Monoethyl Mercury in Aqueous Samples

> Annie Carter and Mi Sun Um, Brooks Rand Labs Seattle Washington USA

> > 0.6

0.5

0.4

(ng/L)

tion



Average RSD: 7%

Std Dev of RSD: 6%

Monomethylmercury (MMHg) has been studied extensively since reliable methods for quantification were developed in the late 1980s. It has been well established that MMHg in the environment is often a byproduct of microbial metabolism, bioaccumulates in the food chain, and is one of the most toxic forms of mercury. Less attention has been given to another organic mercury species, monoethylmercury (EtHg).

One of the sources of EtHg is as a metabolite of thiomerisol, a commonly-used preservative and antimicrobial in many vaccines. As awareness of the toxicity of alkylmercury compounds has increased, so has the concern over their presence in humans and the environment. To study EtHg at lowlevel concentrations, as are expected in most environmental samples, a robust and sensitive method for detection and quantification is required. approved methodology for MMHg analysis in waters, except that NaBPr₄ is used as the derivatization reagent. Additionally, a method has been developed for analyzing water samples for MMHg without a distillation step. In this method, sodium tetraethylborate (NaBEt₄) is added directly to a buffered sample just before analysis for determination by "direct ethylation" (Yu, 2006). The goal of this study was to combine those two techniques into a "direct propylation" method, allowing the simultaneous analysis of MMHg and EtHg in aqueous samples.

Comparing Derivatization and Prep Techniques

An analytical technique utilizing sodium tetra(n-propyl)borate $(NaBPr_4)$ has been developed to determine MMHg and EtHg concentrations in prepared biological samples with great success (Gibičar, 2007). This method is similar to the EPA-

For this direct propylation method, an acidpreserved water sample is buffered just prior to zeroheadspace derivatization with NaBPr₄, followed by cold vapor pre-concentration onto a Tenax[™] trap, isothermal gas chromatography (GC) separation, thermal decomposition, and atomic fluorescence spectrometric detection. All analyses were performed using a Brooks Rand Labs MERX[™] system configured for automated MMHg analysis.



Graph 2: The samples were analyzed for EtHg twice on different days. Both results are graphed. Most samples had acceptable duplicate precision with RPDs of < 20%. Only two samples that were above the EtHg RL had poor duplicate precision. The corresponding MMHg results for those samples had excellent duplicate precision indicating the direct propylation technique is not the cause of the variability seen in these samples results.

Matrix Comparison



Graph 1: To verify the efficiency of the derivatization using NaBPr₄, a side-by-side comparison of the results for three different MMHg analysis methods was performed: *Distillation* – results for distillation followed by derivatization with NaBEt₄, *Direct Ethylation* – results for direct ethylation with NaBEt₄, and *Direct Propylation* – results for direct propylation with NaBPr₄. The RL for the distillation method (0.050 ng/L) was assumed for all techniques.

Graph 3: When MMHg results from direct propylation are compared to the widely-accepted distillation method, the regression is linear with an R² value of 0.975. This shows that both methods give very similar MMHg results, confirming the effectiveness of NaBPr₄ as a derivatization reagent for alkylmercury species in various aqueous matrices.



Graph 4: Three common environmental water matrix types were tested: Fresh *Water* (preserved to 0.4% with HCI), Water + NOM (reagent water high in dissolved organic carbon, spiked with ~10 mg/L of natural organic matter preserved with HCI), and Salt *Water* (preserved with H_2SO_4). All of the test samples were fortified with MMHg and EtHg, since the naturally occurring levels were very low. They were analyzed in triplicate two months after fortification and preservation. They were analyzed at 1, 5, 10, 20, and 40 mL aliquot sizes to



investigate the influence of the matrix on analyte recovery. For the *Fresh Water* and *Water* + *NOM* matrices, both MMHg and EtHg recoveries were \geq 90% at all analysis volumes. MMHg recoveries in the *Salt Water* matrix were biased low (70-87%), even at the lowest aliquot size, and continued to decrease with increasing analysis volumes. EtHg showed extremely low recoveries at all analysis volumes in *Salt Water*. As a result of the low recoveries, additional *Salt Water* samples were fortified with known amounts of EtHg just previous to derivatization, the results are graphed as *Salt Water Spike*. The recoveries for the *Salt Water Spike* samples were good at the 1 mL and 5 mL aliquot sizes, and decreased sharply at greater analysis volumes, indicating some interference from the *Salt Water* matrix. This matrix may be dealkylating the EtHg, which is an occurrence that has been previously reported (Suda, 1993).



Conclusions

Direct propylation of samples was an excellent derivatization method for simultaneous MMHg and EtHg determination.

Graph 5: MDL study of MMHg and EtHg in *Fresh Water* fortified with 0.030 ng/L of each analyte. Results slightly higher than the fortified value are expected due to the ambient level of the analytes in the *Fresh Water* sample.





Graph 6: MDL study of MMHg and EtHg in reagent water fortified with 0.025 ng/L of each analyte.

- There were very few issues caused by matrix interferences.
- The MMHg concentrations determined by direct propylation compared very well to both direct ethylation and distillation prior to ethylation.
- The EtHg analyses demonstrated excellent precision.
- Low method detection limits (MDL) were obtained for both MMHg and EtHg using direct propylation.

Pros

- simultaneous determination of MMHg and EtHg, and potentially Hg(II)
- low MDLs for both MMHg and EtHg
- quick and easy
- good accuracy and precision, as demonstrated through matrix spike recoveries
- can use commercially available instrumentation -Brooks Rand Lab MERX[™]

Cons

- NaBPr₄ is expensive
- matrix effects noted for salt water matrix; potentially other interferences for matrices other than fresh waters
- broad peaks for EtHg and Hg(II) a GC with temperature ramping capabilities would produce sharper peaks

References

- 1. Abranko, L.; Jokai, Z.; Fodor, P. Investigation of the species-specific degradation behavior of methylmercury and ethylmercury under microwave irradiation. Anal. Bioanal. Chem. 2005, 383, 448-453
- 2. Bloom, N. Determination of pictogram levels of methylmercury by aqueous phase ethylation, flowed by cryogenic gas chromatography with cold vapour atomic fluorescence detection. Can. J. Fish. Aquat. Sci. 1989, 46, 1131-11
- 3. De Smaele, T.; Moens, L.; Dams, R.; Sandra, P.; Van der Eycken, J.; Vandyck, J. Sodium tetra(n-propyl)borate: a novel aqueous in situ derivatization reagent for the simultaneous determination of organomercury, -lead and -tin compounds with capillary gas chromatography-inductively coupled plasma mass spectrometry. J. Chromatogr. A. **1998**, 793, 99-106. 4. EPA Draft Method 1630, 2001, Methyl Mercury in Water by Distillation, Aqueous Ethylation, Purge and Trap, and CVAFS.
- 5. Gibičar, D.; Logar, M.; Horvat, N.; Marn-Pernat, A.; Ponkvar, R.; Horvat, M. Simultaneous determination of trace levels of ethylmercury and methylmercury in biological samples and vaccines using sodium tetra(n-propyl)borate as derivatizing agent. Anal. Bioanal. Chem. 2007, 388:329-340.
- 6. Grinberg, P.; Campos, R. C.; Mester, Z.; Sturgeon, R. E. A comparison of alkyl derivitization methods for speciation of mercury based on sold phase microextraction gas chromatography with furnace atomization plasma emission spectrometry detection. J. Anal. At. Spectrom. 2003, 18:902-909.
- 7. Suda, I.; Suda, M.; Hirayama, K.; Degradation of methyl and ethyl mercury by singlet oxygen generated from sea water exposed to sunlight or ultraviolet light. Arch. Toxicol. 1993, 67, 365-368
- 8. Yu, X.; Booeshaghi, F.; Chandraskhar, T. M.; Luo, Y. L.; Tate, K. Determination of Methyl Mercury at Sup-Picogram Levels by Direct Ethylation and an Automated Technique using P&T/GC/AF; Florida Department of Environmental Protection, Bureau of Laboratories, 2006.