

Determination of Methylmercury in Human Hair: A Comparison of Digestion Methods



Annie Carter, Mi Sun Um, and Brittany Nelson
Brooks Rand Labs, Seattle WA USA

Introduction

Mercury is a well-known human health hazard and is routinely monitored in public health studies in order to assess its effects. Total mercury is frequently measured to assess human exposure, due to the relative ease of analysis; however, this information does not always provide sufficient information to accurately assess the source of exposure and/or risk to a person or community.

Because there is no definitive correlation factor establishing the ratio of methylmercury to total mercury in biological samples, it is important to measure methylmercury in human biomonitoring studies. Minimally invasive sample collection is also desirable in most instances; therefore, human blood, urine, and hair are typically collected for analysis. Urine is not an acceptable biomonitoring matrix for methylmercury since methylmercury is nearly completely absorbed through the gut (Gochfeld, 2003). Therefore hair is particularly valuable, as it provides an exposure timeline, rather than just the instantaneous snapshot of body burden provided by blood analysis.

One commonly employed digestion technique for the determination of methylmercury in tissue samples is a potassium hydroxide and methanol (KOH/MeOH) digestion (Bloom 1992). The KOH/MeOH digestion when used for hair is often subject to matrix interferences as noticed by low matrix spike recoveries at analysis and increased results when analyzed at a dilution. Though the samples can typically be diluted 100 fold and exhibit acceptable matrix spike recoveries, this is not a feasible approach as it also increases the method reporting limit (MRL) by a factor of 100.

When the MRL is elevated that high, many of the results are then near or below the MRL, making them estimates and therefore not a useful biomonitoring tool. In this study, an alternative digestion technique using nitric acid (HNO₃) originally developed for the digestion of insects was applied to the digestion of hair samples and then compared to the commonly used KOH/MeOH digestion procedure for the analysis of methylmercury in hair samples (Hammerschmidt, 2005).


Materials and Methods

Analysis was done following EPA Method 1630. The samples were adjusted to a pH of 4.5 - 5, and then ethylated using 0.050 mL of a 1% sodium tetraethylborate in 2% potassium hydroxide solution (Bloom, 1989). The ethylated samples were brought up to zero headspace with reagent water in a septa topped vial and placed on the autosampler for analysis. All analyses were performed using a Brooks Rand Labs MERX methylmercury autoanalyzer.


Sample Preparation & Digestion Methods		
	Potassium Hydroxide in Methanol (KOH/MeOH)	Nitric Acid (HNO ₃)
Sample Mass	100 mg	100 mg
Digestion Solution	25% KOH in Methanol (1 mL)	4M Nitric Acid (20 mL)
Digestion Temperature/Time	65°C for 4 hours	65°C for 24 hours
Diluent	Methanol	None
Final Volume	2.5 mL	20 mL
Analysis Volume	0.030 mL	0.500 mL
MDL/MRL	1.0/3.0 ng/g	0.5/1.5 ng/g

Figure 1: The two digestion schemes used are outlined in this table. The notable differences between the two methods are the digestion solution used and the final volume.

Results



Other Biota Sample Types



Discussion

Two digestion methods were compared: 25% potassium hydroxide in methanol and 4M nitric acid. Both methods are based on procedures described in the peer-reviewed literature and were assessed for accuracy, reliability, and ease of preparation for hair samples. After digestion, all samples were analyzed following EPA Method 1630 (CV-GC-AFS).

The KOH/MeOH digestion method exhibited significant matrix interference, as demonstrated by the recoveries of matrix spikes. Significant dilution of the KOH/MeOH digestion did yield results that were no longer impacted by interference, but the large dilution resulted in a significant increase of the MRL to a level near or above the level of many of the samples.

The nitric acid digestion method did not exhibit a substantial level of matrix interference, and excellent recoveries were achieved for matrix spikes and reference materials even when no dilution was made.

Though the two digestion methods are equally easy to prepare, the HNO₃ method was a more accurate method for the analysis of hair, fur, and feather samples as they were all able to be analyzed at default volume with no significant matrix interferences.

Literature Cited

Bloom, N. S. On the chemical form of mercury in edible fish and marine invertebrate tissue. *Can. J. Fish. Aquat. Sci.* 1992, 49, 1010-1017.

Bloom, N.S. Determination of picogram levels of methylmercury by aqueous phase ethylation, followed by cryogenic gas chromatography with cold vapor atomic fluorescence detection. *Can. J. Fish Aquat. Sci.* 1989, 46, 1131-1140.

EPA Method 1630. Methyl Mercury in Water by Distillation, Aqueous Ethylation, Purge and Trap, and Cold Vapor Atomic Fluorescence Spectrometry, 2001.

Gochfeld, M. Cases of Mercury Exposure, Bioavailability, and Absorption. *Ecotoxicology and Environmental Safety*, 2003, 56, 174-179.

Hammerschmidt, C. R.; Fitzgerald, W. F. Methylmercury in Mosquitoes Related to Atmospheric Mercury Deposition and Contamination. *Environ. Sci. Technol.* 2005, 39, 3034-3039.