



Introduction

Contamination from trace metals, such as lead (Pb), in soils can lead to potential health risks from ingestion of the soil or vegetables grown in the soil. Typical risk assessment models apply an assumed number ranging from 60 - 100% of contaminant metals in a soil sample as bioavailable. These assumptions can result in incorrect health risk assessments or unnecessarily high remediation costs because actual bioavailability can range from 0 - 100%. Currently, in vivo animal studies are often used to determine bioavailability, which are expensive and time consuming. In vivo methods are giving way to new in vitro methods, such as EPA Method 1340 "In Vitro Bioaccessibility Assay for Lead in Soil". The California Department of Toxic Substances Control recently published recommended methodology for evaluating site-specific arsenic (As) bioaccessibility in soils. This method, the California Arsenic Bioaccessibility (CAB) Method, was developed to improve the accuracy of risk assessments and assist regulators in making responsible remediation decisions while still protecting human health.

Bioavailability refers to the concentration of metal that is readily absorbed when ingested and studied in vivo, while bioaccessibility refers to metals extracted by in vitro methods. In vitro bioaccessibility (IVBA) can be defined as the bioaccessible portion of total metal concentration in percentage, see equation below.

IVBA = (concentration of bioaccessible)/(total concentration)

Relative bioavailability (RBA) of Pb and As by in vivo testing has previously been shown to correlate to IVBA results by EPA 1340 and CAB, respectively. The ratio of RBA to IVBA can be defined by the following equations. Each equation is both metal and extraction specific. RBA Pb is calculated using the IVBA Pb recovery from EPA 1340, and RBA As is calculated using IVBA As recovery from CAB.^{1, 2}

> Pb: $RBA = (0.878 \times IVBA) - 0.028$ As: $RBA = (0.81 \times IVBA) + 3.2$

Note: both methods have an upper concentration limit. EPA 1340 is not recommended for total Pb concentrations >50,000 mg/kg, and CAB is not recommended for total As concentrations >1500 mg/kg.

It is well known that different molecular forms of Pb and As have different toxicological effects; however, the correlation between elemental species in the soil and the bioaccessibility of the metals from that soil has not been well established. For instance, inorganic As was used in pesticides and herbicides in US agriculture until 1993, and organic As was used up until 2013, while monosodium methanearsonate (MSMA), is still approved for use in the US as a herbicide in cotton agriculture.³ All of these compounds could be present in contaminated soil and each has a different degree of bioaccessibility.

For this study, bioaccessible As and Pb concentrations in soil samples were determined using both EPA Method 1340 and the CAB Method, respectively, with sample analysis by inductively coupled plasma mass spectrometry utilizing collision and reaction cell technology (ICP-QQQ-MS). Both As and Pb were quantified in both extractions to evaluate if a streamlined method could be developed and applied to both elements. In addition, samples were analyzed for total recoverable As and Pb using EPA Method 3050B to assess the RBA & IVBA of each metal in the samples. In order to evaluate if there was a correlation between the bioaccessible metals concentrations and specific As compounds, direct quantification of different molecular forms of the metal (e.g., arsenite, arsenate, etc.) were performed on the soil samples.

Methods

Samples were extracted for As and Pb by CAB and EPA 1340, then extracts were analyzed by ICP-QQQ-MS on an Agilent 8800. Both extraction methods are designed to simulate human digestion. Soil samples were mixed with an extraction fluid, designed to simulate stomach acid, rotated and incubated at body temperature. Organic compounds in the extraction fluid vary based on the extraction method.

Soil samples were digested for total As and Pb by EPA 3050B, and analyzed by ICP-QQQ-MS on an Agilent 8800. EPA 3050B is a digestion of solid matrices for environmentally available metals. While this digestion will dissolve almost all elements that could become "environmentally available" it is not a true total digestion. Strongly bound analytes, for example silicate bound structures, are not solubilized by this method.

Arsenic species were extracted from the soil samples using two separate (acidic and basic) extractions. With this method, arsenite (As(III)), arsenate (As(V)), monomethyl arsenic acid (MMA), dimethyl arsenic acid (DMA) were directly measured and quantified utilizing ion chromatography inductively coupled plasma mass spectrometry (IC-ICP-MS) analysis on an Agilent 7700. Data from both extracts were evaluated to determine the most appropriate extraction conditions based on the chemistry inherent within each individual sample. Comparison of CAB and EPA 1340 is shown in Figure 1. CAB and EPA 1340 extractions were also performed on a known reference material, Montana Soil II NIST 2711a. Recoveries of the reference material were compared to known IVBA values to confirm proper extraction of the target metal.

Comparison of CAB and EPA 1340 Extraction Methods					
	EPA 1340	САВ			
Initial Sample Mass	1.0 grams	1.0 grams			
Extraction Solution Volume	100 mL	150 mL			
Incubation Temp	37 ± 2 °C	37 ± 2 °C			
Rotation Speed	30 ± 2 rpm	30 ± 2 rpm			
Incubation Time	60 minutes	120 minutes			
Extraction Time	60-90 minutes	120-145 minutes			
Extraction Fluid	Glycine Hydrochloric Acid	Pepsin NaCl Ascorbic Acid Hydrochloric Acid			
Extraction pH	1.5 ± 0.5 pH units	1.50 ± 0.01 pH units			
Number of pH adjustments	1	2			
Number of pH checks	2	3			

Figure 1. Method comparison of CAB and EPA 1340 Extractions. Difference I methods are in bold.

EPA 1340
Dissolve 60.6 g of gl into 1.5 L of DI H
Heat solution to 37 =
pH adjust heated extraction with HCl to 1.5 ± 0.5 p
Bring to final volume of 2.0
Add 100 mL of extraction sol of each dried and sieved s
Incubate samples at 37 ± rotating at 30 ± 2 rpm fo
Remove from incubat allow sediment to s
Filter solution through I
Check pH of each final each should be 1.5 ± 0.5 p
within range analyze / If by ICP-QQQ-MS / re-

Comparison of EPA 1340 and CAB Method for **Bioaccessible Lead and Arsenic in Soils**

Brian Smith, Ariana Dionisio, Annie Carter Brooks Applied Labs, 18804 North Creek Parkway, Suite 100, Bothell, WA 98011 USA





Results

Pb Recovery EPA 1340 vs CAB



Figure 2. Pb extraction from soil samples was performed by EPA 1340 and CAB. Results show consistently higher concentrations of Pb are extracted when using the CAB method. However, the percent increase of Pb relative bioaccessibility by CAB extraction varied greatly between different soil matrices. While some soil samples only showed a slight increase in extracted Pb, sample #3 increased by 10-fold. Known value of reference material, Montana Soil II NIST2711a, confirms proper extraction of bioaccessible Pb by EPA 1340 and bias high recovery by CAB. *R designates a standard reference material, Montana Soil II NIST 2711a

As Recovery EPA 1340 vs. CAB vs. As(V)



Figure 3. Arsenic extraction from soil samples was performed by EPA 1340 and CAB. Results consistently show significantly lower concentrations of As are extracted when using the EPA 1340 method. Known value of reference material confirms proper extraction of IVBA As by CAB and biased low results by EPA 1340.



As can vary greatly regardless of the total As concentration.

As Speciation

Comparison of Arsenic, Bioaccessible Arsenic, and Arsenic Species (mg/kg)						
Sample	3050B	САВ	As(III)	As(V)	DMA	MMA
R	92.7	74.9	0.159	62.2	U	U
1	604	70.7	U	67.9	U	U
2	318	83.5	0.250	120	U	U
3	653	112	U	112	U	U
4	342	195	U	297	0.170	U
5	300	227	U	243	0.100	U
6	187	130	0.204	95.6	U	U
7	1440	1330	1.59	1030	0.950	U

Figure 5. The dominant As species was As(V) in every soil sample tested. Note: As speciation extraction is not a total recoverable method and the sum of the As species is typically less than the As concentration by EPA 3050B. Arsenic bound to some minerals is not recoverable by As speciation extraction however, a more rigorous digestion of mineral bound As is not possible because it will cause species conversion. *R designates known reference material, Montana Soil II NIST 2711a. ** U value in table designates concentration of analyte was under the detect limits.







Figure 7. There is some evidence to suggest a correlation between As(V) and bioaccessible As however. it is not a strong correlation. R² values of 0.9768 with all points included, and 0.814 with a single outlier removed.

Relative Bioavailability of Pb and As



assumed concentration.

Quality Control

Quality control (QC) samples were prepared and analyzed with all extractions and digests. QC samples were analyzed in the same manner as sample.

Method Blank 1	
Method Blank 2	
Method Blank 3	
Blank Spike	
Reference Material	
Source - Duplicate	
Duplicate	
Source - Matrix Spikes	
Matrix Spike**	
Matrix Spike Duplicate**	
Method Blank 1	
Method Blank 2	
Method Blank 3	
Method Blank 3 Blank Spike	
Method Blank 3 Blank Spike Reference Material	
Method Blank 3 Blank Spike Reference Material Source - Duplicate	
Method Blank 3 Blank Spike Reference Material Source - Duplicate Duplicate	
Method Blank 3 Blank Spike Reference Material Source - Duplicate Duplicate Source - Matrix Spikes	
Method Blank 3 Blank Spike Reference Material Source - Duplicate Duplicate Source - Matrix Spikes Matrix Spike**	

CAB and EPA 1340 confirm samples were properly extracted for target analytes Sample duplicates noted very low relative percent differences (RPD) in all cases, <4%. Blank spikes recovery very close to 100% in all instances. Reference material recoveries confirm EPA 1340 under recovers for As and CAB over recovers for Pb. *U designates results that were under the detection limits of the instrument.** Matrix spikes (MS) and matrix spike duplicates (MSD) were spiked at the same level as blank spikes. This resulted in under-spiked MS and MSD, poor spiking recoveries are expected.



The Ohio State University

Figure 8. Relative Bioavailability (RBA) was calculated from IVBA for both As and Pb. RBA As was calculated using the CAB extraction recoveries, and RBA Pb was calculated using the 1340 extraction recoveries. Assumed RBA typically ranges from 60-100%. In most cases RBA concentrations are much lower than

Quality Control						
CAB Extraction Quality Control Sample Recoveries						
As (mg/kg)	Percent Recovery	Relative Percent Difference	Pb (mg/kg)	Percent Recovery	Relative Percent Difference	
0.215			0.213			
0.214			0.111			
0.181			0.114			
10.8	108%		1.06	106%		
74.9			1290	116%		
227			1910			
234		3.00%	1850		3.00%	
195			1580			
200	46.0%**		1600	2190%**		
203	78.0%**	52.0%**	1610	2740%**	61.0%**	

340 Extraction Quality Control Sample Recoveries

As (mg/kg)	Percent Recovery	Relative Percent Difference	Pb (mg/kg)	Percent Recovery	Relative Percent Difference
0.002			0.01		
0.0004			0.004		
U*			0.007		
9.96	100%		1.04	104%	
54.1	70.0%		1220	111%	
123			1820		
121		1.00%	1760		3.00%
467			3230		
475	79.0%**		3170	-6230%**	
472	57.0%**	33.0%**	3200	-3200%**	64.0%**

Discussion

Two extraction methods for bioaccessible As and Pb were performed on various soil samples: EPA 1340 for bioaccessible Pb, and CAB for bioaccessible As. Samples were also analyzed by EPA 3050B for total recoverable metals and for various As species by IC-ICP-MS.

The CAB method showed higher concentrations of both As and Pb than EPA Method 1340. This indicates a higher extraction efficiency for the CAB method. This is likely due to the different reagents used in the extraction fluids and the increased incubation and extraction times in the CAB method. As a result, using the CAB method for Pb would provide results that are biased high for bioaccessible Pb for many samples. Inversely, using the EPA 1340 method for As would produce results that are biased low for bioaccessible As.

In all analyzed soil samples As(V) made up the vast majority of the detected As species. Some correlation was seen between IVBA As and As(V) however, relative percent differences ranged from 0-42%. IVBA As and As(V) correlation may be depend on each specific soil matrix.

Conclusion

• Varying extracted concentrations in each method illustrate the importance of using the appropriate extraction method for the desired target analyte.

 The newly developed CAB method gave consistent results for bioaccessible As but yields concentrations for bioaccessible Pb that are biased high compared to EPA 1340

• EPA 1340 yielded expected results for bioaccessible Pb but bias low results for bioaccessible As compared to CAB.

• As(V) speciation may give some indication of bioaccessible As but should not be used as a proxy for the CAB extraction method.

 Brooks Applied Labs, LLC plans to investigate applying these bioaccessibility extractions to other metals such as cadmium and selenium, due to their similar properties to lead and arsenic.

Acknowledgements

The authors of this poster would like to thank Shane Whitacre and Nicholas Basta of The Ohio State University (OSU) and Valerie Hanley of The California Department of Toxic Substance Control (DTSC) for their collaboration on this project. Additional thanks to DTSC and OSU for hosting the round robin study of The California Arsenic Bioaccessibility (CAB) Method.

References

1. Validated Test Method 1340: In Vitro Bioaccessibility Assay for Lead in Soil. (2017, March 31). Retrieved June 1, 2017, from https://www.epa.gov/hw-sw846/validated-test-method-1340-<u>vitro-bioaccessibility-assay-lead-soil</u>

2. Whitacre, S., Basta, N., Stevens, B., Hanley, V., Anderson, R., & Scheckel, K. (2017). Modification of an existing in vitro method to predict relative bioavailable arsenic in soils. Chemosphere, 180, 545-552. doi:10.1016/j.chemosphere.2017.03.134

3. IARC Working Group on the Evaluation of Carcinogenic Risk to Humans. Arsenic, Metals, Fibres and Dusts. Lyon (FR): International Agency for Research on Cancer; 2012. (IARC Monographs on the Evaluation of Carcinogenic Risks to Humans, No. 100C.) ARSENIC AND ARSENIC COMPOUNDS. Available from:

<u>https://www.ncbi.nlm.nih.gov/books/NBK304380/</u>